



GH3+ Scale module (accessories)

The scale module is very useful anywhere the precise weighing of an individual is needed to monitor his or her changes in weight.

Scale module

- The scale module can be used for weighing patients whose weight needs to be monitored.
- Medical scales can and may only be calibrated by an accredited individual/company, and Guldmann's integrated scales are therefore not medical scales.
- Guldmann's scale module can be calibrated by a technician approved by Guldmann (local technical staff).
- The accuracy of Guldmann's integrated scale is the same as medically certified scales.
- The GH3+ solution complies with the EN/ISO 10535:2007 standard.

Specifications of the scale module

Capacity	0 - SWL
Accuracy	+/- 0.1 % at max load
Display resolution (d)	0.1 kg (0.22 lbs)
Accuracy of repetition	< 0.1 kg at 0-200 kg (0-440 lbs)
	< 0.2 kg at 0-350 kg (0-770 lbs)
Minimum weight	5 kg (11 lbs)
Display type	LCD in Hand Control